



DROPPING FUNNELS, CYLINDRICAL, WITHOUT PRESSURE BALANCE

Borosilicate 3.3 glass.

- Ground glass cone for flasks and assemblies
- PTFE key stopcock and PP stopper
- Very good chemical resistance
- High temperature resistance
- Minimal thermal expansion giving relatively high resistance to temperature changes

PTFE key stopcock for non stick, grease-free performance.

ISO 4800

Capacity (ml)	Cone size (NS)	Socket size (NS)	Pk	Cat. No.
50	14/23	14/23	1	201-1238
50	29/32	29/32	1	201-1239
100	14/23	14/23	1	201-1240
100	29/32	29/32	1	201-1241
250	29/32	29/32	1	201-1242
500	29/32	29/32	1	201-1243
1000	29/32	29/32	1	201-1244



DROPPING FUNNELS, CYLINDRICAL, WITH PRESSURE BALANCE

Borosilicate 3.3 glass.

- Standard tapered joints top and bottom, drip stem and pressure-equalising arm
- PTFE key stopcock and plastic stopper
- Very good chemical resistance
- Very good temperature resistance

Drip stem does not extend beyond joint, thereby eliminating a possible cause of breakage. PTFE key for grease-free performance.

Capacity (ml)	Cone size (NS)	Socket size (NS)	Pk	Cat. No.
50	14/23	14/23	1	201-1245
50	29/32	29/32	1	201-1246
100	14/23	14/23	1	201-1247
100	29/32	29/32	1	201-1248
250	29/32	29/32	1	201-1249
500	29/32	29/32	1	201-1250



FILTER HOLDER

Glass filter holder with stainless steel screen, 47 mm.

Description	Pk	Cat. No.
Glass filter holder with stainless steel screen	1	511-0666



FILTRATION APPARATUS

Borosilicate 3.3 glass.

- Very good chemical resistance
- Very good temperature resistance
- Available with 1000 or 2000 ml flask
- Funnel with 300 or 500 ml capacity

All-glass filtration apparatus, with funnel, fritted base, cap, 47 mm Ø spring clamp and ground joint flask. This unit is recommended for general filtration. Applications include filtering corrosive solutions, particulate removal from HPLC solvents and filtration analysis with solid phase extraction discs. Support base has a sealed-in, coarse porosity, fritted disc.

Capacity (ml)	Description	Pk	Cat. No.
300	With 1000 ml flask	1	511-0265
500	With 2000 ml flask	1	511-0266

Description	Pk	Cat. No.
Accessories		
Spare glass base and cap for filtration apparatus, 47 mm	1	511-0270
Spare glass funnel 300 ml for filtration apparatus, borosilicate 3.3 glass	1	511-0267
Spare glass funnel 500 ml for filtration apparatus, borosilicate 3.3 glass	1	511-0268
Spare ground joint flask for filtration apparatus, 1000 ml capacity	1	511-0271
Spare ground joint flask for filtration apparatus, 2000 ml capacity	1	511-0272
Spring clamp, aluminium, 47 mm for filtration apparatus	1	511-0269



FILTERING FLASKS, WITH GLASS HOSE CONNECTION

Borosilicate 3.3 glass.

- Suitable for hot and cold filtration
- Highly resistant to chemicals and thermal shock
- High temperature resistance
- White labelling area

Due to the heavy wall thickness the flask is vacuum tight.

Capacity (ml)	Pk	Cat. No.
100	1	511-0260
250	1	511-0261
500	1	511-0262
1000	1	511-0263
2000	1	511-0264



Solutions for every day use

VWR for weighing

- Balances - top pan, analytical and micro
- Weighing boats and scoops
- Calibration - servicing, maintenance and everyday calibration weights

For your copy
visit vwr.com