

- Equipment
- ErgoOne® Pipettes
- Consumables
- **Sample Storage**



# Sample Storage.

# Sample storage solutions.



## TwisStar Rack

Four, different coloured polypropylene sections that can be locked together to suit your needs. Holds an array of different tubes and a PCR plate, it's ideal for use in the hood! Temperature use: -90°C to 121°C.

Art. No. E2344-1099

**OFFER:**



## PCR Lock-Rack

Polypropylene racks which lock together to make larger racks. Transparent lid fits over two sections locked together. Each section holds 24 x 1.5 / 2.0 ml tubes, 16 x 0.5 ml tubes and 60 x 0.2 ml tubes. Temperature use: -90°C to 121°C.

Art. No. E2344-2099

**OFFER:**



## Polycarbonate Cryogenic Storage Boxes


White base with transparent, lift-off lid with printed reference grid. Temperature use from -196°C to 121°C. Standard 130x130 footprint for use in standard freezer racking.

Art. No. I2381-2305 (81-place)

**OFFER:**

Art. No. I2510-5005 (100-place)

**OFFER:**




## StarChill MCT, Green/Yellow

24-well rack with lid. Use with 0.5, 1.5 and 2.0 ml microcentrifuge / screw cap tubes. Temperature colour change at 7° C.

Art. No. E2396-5636

**OFFER:**




## StarChill PCR, Purple/Pink

96-well rack with lid. Use with 0.2 ml tubes or strips, or a 0.2 ml 96-well PCR plate. Temperature colour change at 7° C.

Art. No. E2396-2004

**OFFER:**



## StarChill SBS

96-well rack (no lid), for use on AUTOMATED SYSTEMS with 0.2 ml tubes or strips, or a 0.2 ml 96-well PCR plate.

Art. No. E2396-5656

**OFFER:**